

Assessment of Mitochondrial Respiration During Hypothermic Storage of Liver Biopsies Following Normothermic Machine Perfusion

Julia Hofmann^{1†}, Alexander Kofler^{1†}, Melanie Schartner¹, Madita L. Buch¹, Martin Hermann¹, Bettina Zelger², Dietmar Öfner¹, Rupert Oberhuber¹, Theresa Hautz¹, Stefan Schneeberger¹ and Andras T. Meszaros^{1*}

¹OrganLife Laboratory, Department of Visceral, Transplant and Thoracic Surgery, Medical University of Innsbruck, Innsbruck, Austria, ²Institute of Pathology, Neuropathology and Molecular Pathology, Medical University of Innsbruck, Innsbruck, Austria

Organ quality can be assessed prior to transplantation, during normothermic machine perfusion (NMP) of the liver. Evaluation of mitochondrial function by high-resolution respirometry (HRR) may serve as a viability assessment concept in this setting. Freshly collected tissue is considered as optimal sample for HRR, but due to technical and personnel requirements, more flexible and schedulable measurements are needed. However, the impact of cold storage following NMP before processing biopsy samples for mitochondrial analysis remains unknown. We aimed at establishing an appropriate storage protocol of liver biopsies for HRR. Wedge biopsies of 5 human livers during NMP were obtained and assessed by HRR. Analysis was performed after 0, 4, 8, and 12 h of hypothermic storage (HTS) in HTK organ preservation solution at 4°C. With HTS up to 4 h, mitochondrial performance did not decrease in HTS samples compared with 0 h (OXPHOS, 44.62 [34.75-60.15] pmol·s⁻¹·mg wet mass⁻¹ vs. 43.73 [40.69-57.71], median [IQR], p > 0.999). However, at HTS beyond 4 h, mitochondrial respiration decreased. We conclude that HTS can be safely applied for extending the biopsy measurement window for up to 4 h to determine organ quality, but also that human liver respiration degrades beyond 4 h HTS following NMP.

OPEN ACCESS

*Correspondence

ransplant

ternational

Andras T. Meszaros, ⊠ andras.meszaros@i-med.ac.at

[†]These authors have contributed equally to this work and share first authorship

> Received: 31 January 2024 Accepted: 07 May 2024 Published: 23 May 2024

Citation:

Hofmann J, Kofler A, Schartner M, Buch ML, Hermann M, Zelger B, Öfner D, Oberhuber R, Hautz T, Schneeberger S and Meszaros AT (2024) Assessment of Mitochondrial Respiration During Hypothermic Storage of Liver Biopsies Following Normothermic Machine Perfusion. Transpl Int 37:12787. doi: 10.3389/ti.2024.12787 Keywords: liver transplantation, mitochondria, normothermic machine perfusion, high-resolution respirometry, hypothermic storage

INTRODUCTION

The serious shortage of donor organs is a daily challenge in liver transplantation (LT). To overcome this limitation and reduce waiting list mortality, extended criteria donor (ECD) organs are increasingly utilized for LT [1]. However, a substantial proportion of retrieved organs are not transplanted since they are considered inferior in their macroscopic appearance and expected function. Thus, it is essential to evaluate the quality of these grafts and predict the clinical outcome with the help of reliable biomarkers [2]. Machine perfusion is an emerging organ preservation strategy. Donor livers are perfused under normothermic conditions, which also allows for assessment of the organ function prior to transplantation. For this purpose, several biomarkers have been recently suggested. Transplant centers usually rely on a combination of different parameters for decision-making, including lactate clearance, level of perfusate transaminases, bile production and its



composition. However, since robustness of those parameters with regard to their predictive capacity is limited, new biomarkers are under investigation [3–5]. Recently, our group demonstrated that mitochondrial bioenergetics assessed by high-resolution respirometry (HRR) in liver biopsies during the first 6 h of NMP predicts graft function after LT [6]. Different methods to analyze mitochondrial status are available, including monitoring free flavin mononucleotide [7–10], a mitochondrial damage parameter. In contrast, HRR allows for evaluation of the actual bioenergetic capacity by assessing mitochondrial respiration. Oxygen consumption of mitochondria is modulated by various substrates and inhibitors and is measured in real-time, thus enabling a comprehensive characterization of energy production, mitochondrial (dys) function and metabolic pathways.

In the above study, we performed HRR in tissue homogenates of freshly collected biopsies. However, this might be challenging in the clinical routine, especially in serial biopsies requiring immediate measurements, each taking approximately 2 h. To reduce this logistical and technical burden, extended hypothermic storage (HTS) would require little effort and ease the logistics. HRR analyses of tissue samples collected in the early phase of NMP could be performed in parallel, and the results would be available in a relevant period of the decision-making process for transplantation [11, 12]. The effect of cold ischemia on mitochondrial respiration has been extensively studied in endothelial cells [13], in heart [14], and in small samples of porcine livers [15]. However, no data on preservation protocols for liver biopsies taken during NMP exist. In this study we aimed to investigate the effect of HTS on mitochondrial respiration and coupling control in clinical liver NMP.

MATERIAL AND METHODS

Study Design and Ethics

Tissue samples of liver allografts (N = 5) subjected to NMP (Metra[®], OrganOx Limited, Oxford, United Kingdom) prior to transplantation were included in this single center experimental study between January and August 2023. The study was approved by the institutional ethics committee (Ethics Commission of the Medical University of Innsbruck, EK Nr. 1175/2018). NMP was performed in accordance with our center's protocols as described by Cardini et al. [16, 17]. After organ acceptance and indication for NMP subjection, livers were included in the study on a rolling-basis when (i) logistics could be arranged and repetitive high-resolution respirometry (ii) measurements during the next 12 h could be ensured. At the time point of inclusion, the performance of the livers during NMP was still under investigation. Wedge biopsies $(20 \text{ mg} \pm 2 \text{ mg wet mass})$ were taken during clinical NMP and divided into 4 approximately equal parts (~5 mg each). The biopsy samples were then directly placed in ice-cold HTK solution (Custodiol[®], Dr. Franz Köhler Chemie GmbH, Bensheim, Germany) and transferred to the laboratory at



FIGURE 1 Representative trace of a high-resolution respirometry measurement applying the SUIT (i) protocol. After addition of the tissue homogenate (0.5 mL), residual oxygen consumption (ROX) was assessed in the absence of fuel substrates and adenosine diphosphate (ADP). Next, ADP (5 μ M) was added, followed by malate (0.1 mM) and octanoylcarnitine (0.5 mM), which allows for assessment of the F_P OXPHOS capacity. Integrity of the outer mitochondrial membrane was tested by titration of cytochrome *c* (10 μ M). The subsequent titration of N_P-linked substrates (2 mM malate, 5 mM pyruvate and 10 mM glutamate) enables evaluation of combined F_P- and N_P -linked OXPHOS capacity was measured by the addition of succinate (10 mM), followed by the titration of Complex I inhibitor rotenone (0.5 μ M), which allows the measurement of S_P-linked OXPHOS capacity only. Lastly, the Complex III inhibitor antimycin A (2.5 μ M) was added, which allows for evaluation of ROX. OXPHOS, oxidative phosphorylation; F_P, fatty acid oxidation-linked pathway; N_P, nicotinamide adenine dinucleotide-linked pathway; S_P, succinate-linked pathway; ROX, residual oxygen consumption.

4°C on ice for further processing. The first biopsy was immediately analyzed for baseline mitochondrial respiratory function values. The further samples were stored at 4°C on ice and analyzed after 4, 8, or 12 h of HTS.

High-Resolution Respirometry (HRR)

Mitochondrial respiration was assessed by HRR (O2k, Oroboros Instruments, Innsbruck, Austria) with devices equipped with two identical 0.5 mL measurement chambers at 37°C and under continuous stirring at 750 rpm [18]. Measurements were preceded by air-calibration of the devices using the mitochondrial respiration medium MiR05 (MiR05-Kit, Oroboros Instruments, Innsbruck, Austria) containing 0.5 mM EGTA, 3 mM MgCl₂ • 6 H₂O, 60 mM lactobionic acid, 20 mM taurine, 10 mM KH₂PO₄, 20 mM HEPES, 110 mM D-sucrose, $1 \text{ g} \cdot \text{L}^{-1}$ essentially fatty acid free bovine serum albumin). Liver wedge biopsies were weighted in a 4°C environment, and subsequently homogenized in 4°C MiR05 using a PBI-Shredder O2k-Set (Oroboros Instruments, Innsbruck, Austria), to disrupt cellular structures and extract mitochondria. The crude tissue homogenate, containing freely accessible mitochondria, was adjusted to a final concentration of $1 \text{ mg wet mass} \cdot \text{mL}^{-1}$, and 500 µL was transferred into the O2k chambers. Reagents were added using Hamilton glass microsyringes (Hamilton Bonaduz AG, Bonaduz, Switzerland) following pre-defined Substrate-Uncoupler-Inhibitor Titration (SUIT) protocols. All measurements were performed in technical duplicates [6]. Oxygen concentration and flux were recorded at an interval of 2 s using DatLab 7 software (Datlab 7.4, Oroboros Instruments, Innsbruck, Austria).

Substrate-Uncoupler-Inhibitor Titration (SUIT) Protocols

We applied two different SUIT protocols. SUIT(i), corresponding SUIT-025 in the online protocol library Bioblast¹, (Figure 1) allows for evaluation of the substrate pathway control in the OXPHOS state. After addition of ADP in saturating concentrations (5 mM), the substrates for the fatty acid oxidation-linked pathway (F_p) , malate (0.1 mM)and octanoylcarnitine (0.5 mM) were titrated, followed by cytochrome c (10 μ M) for the evaluation of the integrity of the mitochondrial outer membrane. Next, substrates for the NADHlinked pathway (N_P) ; pyruvate (5 mM), malate (2 mM) and glutamate (10 mM) were added. The overall capacity of convergent electron transfer pathways was determined after succinate addition (10 mM) (FNS_p). For determination of the succinate-linked pathway alone, (S_p) rotenone $(0.5 \,\mu\text{M})$ was titrated. Finally, to inhibit Complex III and determine residual oxygen consumption (ROX), antimycin A (2.5 µM) was added.

SUIT(ii), corresponding SUIT-006 in the online protocol library Bioblast² (**Figure 2**) provides detailed information on the S_{*P*}-linked mitochondrial respiration. Rotenone (0.5 μ M) was added for inhibition of the mitochondrial Complex I, which was followed by succinate (10 mM) titration. Next, ADP (5 mM) for evaluation of the OXPHOS capacity was added followed by cytochrome *c* (10 μ M) for the evaluation of the integrity of the mitochondrial outer membrane. Carbonyl cyanide

¹https://wiki.oroboros.at/index.php/SUIT-025

²https://wiki.oroboros.at/index.php/SUIT-006



and after the subsequent addition of succinate (10 mM), LEAK respiration could be measured. The addition of ADP (5 mM) enables S_P-linked OXPHOS capacity assessment and with the subsequent titration of cytochrome *c* (10 µM) the integrity of the outer mitochondrial membrane can be evaluated. To assess the electron transfer (ET) capacity, step-wise titration (0.5 µM per step) of CCCP until maximum respiration was performed. Lastly, the Complex III inhibitor antimycin A (2.5 µM) was added, inducing ROX. OXPHOS, oxidative phosphorylation; S_P, succinate-linked pathway; CCCP, carbonyl cyanide p-trifluoromethoxyphenyl hydrazone.

p-trifluoromethoxyphenyl hydrazone (CCCP) was titrated stepwise (0.5 μ M per step) to assess non-coupled respiration. Finally, antimycin A (2.5 μ M) was added to achieve ROX.

Mitochondrial Coupling Control States

Respiration rates were expressed as O_2 flux per wet mass tissue [pmol· $O_2 \cdot s^{-1} \cdot mg^{-1}$]. In addition, respiratory capacities were normalized to an internal reference state for each measurement to determine flux control ratios (*FCR*) for evaluation of SUIT(i). For SUIT(ii), the coupling states LEAK (*L*), OXPHOS (*P*), OXPHOS(*c*) (*P_c*), and Electron Transfer capacity (*ET*) were evaluated. Based on the above, the following control efficiencies were calculated: *P*-*L* control efficiency (1-*L*/*P*), to evaluate the efficiency of the ATP production within S_{*P*}, and cytochrome *c* control efficiency (1-*P*/*P_c*) to evaluate the damage to the mitochondrial outer membrane [17, 19].

Real-Time Confocal Microscopy (RTCM)

Tissue integrity and cell viability of the obtained biopsies were assessed by real-time confocal microscopy (RTCM) as previously described [20, 21]. Briefly, biopsies were incubated with the fluorescent stains Hoechst 33,342 (blue, staining all nuclei, Molecular Probes; final concentration $5 \,\mu \text{g} \cdot \text{mL}^{-1}$), propidium iodide (PI; red, staining the nuclei of dead cells, Molecular Probes; final concentration 500 nM), and wheat germ agglutinin (WGA; green, staining of cell morphology, Molecular Probes, Eugene, OR, United States; $5 \,\mu \text{g} \cdot \text{mL}^{-1}$ final concentration) at room temperature for 15 min. Image acquisition was performed with an UltraVIEW VoX spinning disk confocal system (Perkin Elmer, Wellesley, MA) connected to a Zeiss AxioObserver Z1 microscope (Zeiss, Oberkochen, Germany). Quantification relied on a semi-quantitative scoring system (0-1-2): 0—high number of viable cells, 1—an equal count of viable and dead cells, and 2—low number of viable cells [6].

Histopathology

Liver biopsies were fixed in 10% formalin for 24 h, then embedded in paraffin and cut in 3 μ M sections. Hematoxylin and eosin (H&E) staining was performed following standard protocols to assess necrosis, inflammation, steatosis, fibrosis, and vascular changes. For scoring a modified Suzuki score was applied [6].

Statistical Analysis

Shapiro-Wilk test was employed to examine the normality of the data. Continuous variables are represented as median [interquartile range (IQR)]. Categorial variables are expressed as percentage (%). To assess differences compared to baseline values (0 h sample), Friedman's test with a post-hoc comparison to examine specific pairwise differences by the Dunn's method was performed. Statistical analysis was performed using GraphPad Prism 10 and a *p*-value of <0.05 was considered statistically significant.

RESULTS

Donor and Preservation Characteristics

5 human livers undergoing NMP were included in this study. Four grafts (80%) originated from donors after brain death (DBD). The median age of the donors was 72 [54–78] years, with a median body mass index of 29 [23–33] kg·m⁻². The median donor risk index (DRI) was 2.04 [1.94–2.37], and 80% of organs were procured from ECD. One liver was found suitable

TABLE 1 | Donor and preservation characteristics.

Donor characteristics	Liver 1	Liver 2	Liver 3	Liver 4	Liver 5
Age (years)	65	43	72	79	77
Sex	М	М	М	Μ	М
Body mass index (kg⋅m ⁻²)	29	23	37	29	23
Cause of death	CVA	CVA	CVA	CVA	CVA
Donor type	DBD	DCD	DBD	DBD	DBD
ECD status	no	yes	yes	yes	yes
DRI	1.92	2.43	2.04	2.31	1.96
Steatosis ^a	mild	no	mild	no	no
Fibrosis ^a	mild	no	no	no	no
Hepatitis ^a	minimal	no	no	no	no
Preservation	Liver 1	Liver 2	Liver 3	Liver 4	Liver 5
CIT (min)	538	599	492	306	365
Functional WIT ^b (DCD; min)		18			
Lactate ≤2.5 mmol·L ⁻¹ after 2 h NMP	yes	yes	yes	yes	no
AST peak (NMP) (U·L ⁻¹)	2,417	11,221	1,048	2,159	12,800
ALT peak (NMP) (U·L ⁻¹)	2,588	3,327	909	1,515	11,655
LDH peak (NMP) (U·L ⁻¹)	5,596	16,339	3,036	3,270	25,507
Bile production during NMP	no	yes	yes	yes	yes
Reason for decline	pH instable, no bile	Inhomogeneous perfusion after	Transplanted	Recipient tested positive for	Late lactate clearance,
	production	complex arterial reconstruction		COVID-19, liver was	cholestasis, high LDH
				reallocated	-

^aAssessed in frozen tissue sections.

^b In the context of DCD donors, functional warm ischemia time refers to the interval starting when the mean arterial pressure falls below 50 mmHg or arterial saturation drops below 80%, and ending when cold perfusion starts.

CVA, cerebrovascular accident; ECD, extended criteria donor; DBD, donation after brain death; DCD, donation after circulatory death; DRI, donor risk index; CIT, cold ischemic time; WIT, warm ischemic time; NMP, normothermic machine perfusion; AST, aspartate aminotransferase; ALT, alanine aminotransferase; LDH, lactate dehydrogenase.

for transplantation after NMP, while four livers deemed unsuitable and had to be declined. The individual donor and preservation characteristics, viability parameters and reasons for decline can be found in **Table 1**.

High-Resolution Respirometry for Assessment of Hypothermically Stored Liver Biopsies

With the progression of hypothermic storage, gradual, time-dependent changes in mitochondrial respiration were observed (**Figure 3**).

Characterization of Pathway Contribution

We investigated the NADH- (N_P), FAO- (F_P), and succinatelinked (S_P) oxidative phosphorylation (**Figure 4; Table 2**). Respiration was highest for S_P and lower for N_P and F_P. Concerning S_P, fresh and 4 h HTS measurements revealed similar respiration rates (53.59 [49.51–65.87] pmol·s⁻¹·mg wet mass⁻¹ vs. 51.01 [38.31–66.91] pmol·s⁻¹·mg wet mass⁻¹, p > 0.999), while after 8 h HTS a marked decrease in respiration was observed (43.32 [27.22–55.48] pmol·s⁻¹·mg wet mass⁻¹, p = 0.811). At 12 h, the respiration was even lower: (37.99 [19.61–47.79] pmol·s⁻¹·mg wet mass⁻¹, p = 0.06). In line with this, N_P and F_P remained stable for 4 h of HTS and declined after 8 h: N_P decreased to 3.89 [3.10–4.94] pmol·s⁻¹·mg wet mass⁻¹ after 12 h (p = 0.006). Furthermore, the flux control ratios for each pathway were calculated, to take advantage of the internal normalization relative to the maximal convergent electron pathway respiration (**Figure 4; Table 2**). In the present cohort, in fresh biopsies FCRs of 0.27 and 0.14 were found for F_P and N_P , respectively, which also remained constant over storage time. For S_P a FCR of 1.0 was observed for fresh biopsies and all analyzed HTS biopsies, which indicated that succinate alone is sufficient to saturate OXPHOS capacity.

Succinate-Linked Mitochondrial Respiration

The succinate-linked respiration was analyzed in greater detail as it relates to the prominent substrate pathway in the present and previous studies (Table 3; Figure 5). Measurement of LEAK respiration, representing the non-phosphorylating component of mitochondrial respiration, yielded similar results between fresh tissue biopsies and those measured after 4 h HTS (6.75 [6.69-8.50] $\text{pmol}\cdot\text{s}^{-1}\cdot\text{mg}$ wet mass⁻¹ vs. 7.24 [7.03-8.09] pmol·s⁻¹·mg wet mass⁻¹, p > 0.999). In line with this, OXPHOS capacity, measured with saturating substrate, ADP, and inorganic phosphate concentrations remained stable for 4 h HTS (fresh 43.73 [40.69-57.71] pmol·s⁻¹·mg wet mass⁻¹ vs. 4 h HTS 44.62 [34.75–60.15] pmol·s⁻¹·mg wet mass⁻¹, p > 0.999). The same pattern was also observed for the electron transfer capacity in the noncoupled state (fresh 82.64 [78.22-93.51] pmol·s⁻¹·mg wet mass⁻¹ vs. 4 h HTS 71.09 [68.24–96.07] pmol·s⁻¹·mg wet $mass^{-1}$, p > 0.999). All assessed parameters showed decrease after 8 h HTS and beyond. However, significance was not reached.



FIGURE 3 Assessment of mitochondrial respiration in tissue biopsies after 0, 4, 8, and 12 h of hypothermic storage. (A) ROX was measured after addition of tissue homogenate and was followed by titration of ADP to evaluate OXPHOS capacities for the F, FN and FNS pathway states. (B) ROX was measured after addition of tissue homogenate followed by Complex I inhibition and succinate titration to evaluate S-linked LEAK respiration. ADP and cytochrome *c* addition allowed for measurement of OXPHOS capacity and OXPHOS (c) capacity, respectively. Uncoupler titrations were performed to analyze ET capacity. ADP, adenosine diphosphate; F, fatty acid oxidation; FN, fatty acid oxidation and NADH; FNS, fatty acid oxidation-NADH-succinate; OXPHOS, oxidative phosphorylation; ET, electron transfer; ROX, residual oxygen consumption; thom, tissue homogenate; Oct, octanoylcarnitine; M, malate; P, pyruvate; G, glutamate; S, succinate; Rot, rotenone; Ama, antimycin A; cyt c, cytochrome *c*; CCCP, Carbonyl cyanide p-trifluoromethoxyphenyl hydrazone. Data are represented as median and interguartile range.





TABLE 2 | Pathway-specific oxidative phosphorylation.

Mitochondrial respiration (pmol·s ⁻¹ ·mg wet mass ⁻¹)	0 h	4 h	<i>p</i> -value	8 h	p-value	12 h	<i>p</i> -value
NADH (N _P)	7.84 [6.65–11.78]	6.30 [4.95-8.66]	0.981	5.86 [4.52-7.88]	>0.999	3.89 [3.10-4.94]	0.006
FAO (F _P)	17.20 [11.00-26.18]	14.17 [13.08-27.03]	>0.999	15.40 [9.26–15.89]	0.534	14.99 [6.31–18.38]	0.334
Succinate (S _P)	53.59 [49.51-65.87]	51.01 [38.31–66.91]	>0.999	43.32 [27.22–55.48]	0.811	37.99 [19.61–47.79]	0.06
FCR (relative to FNS)	0 h	4 h	p-value	8 h	p-value	12 h	p-value
NADH (N _P)	0.15 [0.13–0.18]	0.13 [0.10-0.15]	0.662	0.15 [0.12-0.19]	>0.999	0.11 [0.10-0.17]	>0.999
FAO (F _P)	0.26 [0.17-0.35]	0.30 [0.21-0.34]	0.534	0.25 [0.22-0.30]	>0.999	0.27 [0.23-0.35]	0.259
Succinate (SP)	1.00 [0.99–1.00]	0.99 [0.99–0.99]	>0.999	1.03 [1.01–1.08]	0.199	1.03 [1.02–1.07]	0.150

NADH, nicotinamide adenine dinucleotide hydrogen; FAO, fatty acid oxidation; FCR, flux control ratio; FNS, fatty acid oxidation-NADH-succinate. Data are represented as median, interquartile range and min-max. Friedman's test with Dunn's multiple comparisons was applied to assess differences compared to 0 h samples.

TABLE 3 | Mitochondrial respiratory capacities of the S-pathway.

Coupling states (pmol·s ^{-1.} mg wet mass ⁻¹)	0 h	4 h	<i>p</i> -value	8 h	<i>p</i> -value	12 h	<i>p</i> -value
LEAK	6.75 [6.69–8.50]	7.24 [7.03-8.09]	>0.999	6.19 [5.33–8.30]	>0.999	5.51 [5.02-8.54]	>0.999
OXPHOS	43.73 [40.69–57.71	44.62 [34.75-60.15]	>0.999	36.56 [24.38-49.41]	0.811	35.29 [16.18-46.77]	0.199
OXPHOS(c) ET	46.31 [44.55–60.41 82.64 [78.22–93.51]	49.42 [36.75–61.42] 71.09 [68.24–96.07]	>0.999 >0.999	41.17 [29.47–53.69] 60.95 [51.07–86.91]	0.811 0.811	39.71 [21.42–50.43] 52.33 [34.26–75.40]	0.199 0.06

OXPHOS, oxidative phosphorylation; OXPHOS(c), OXPHOS capacity after cytochrome c addition; ET, electron transfer. Data are represented as median, interquartile range and min-max. Friedman's test with Dunn's multiple comparisons was applied to assess differences compared to 0 h samples.

Coupling Control Efficiencies

The calculated mitochondrial coupling control efficiencies of the succinate-pathway did not show significant alterations during HTS up to 12 h compared to fresh biopsies (**Table 4**). Calculation of the *P-L* control efficiency allows for evaluation of the efficiency of ATP production. It was sufficiently high in fresh biopsies (0.85 [0.83–0.86]) but a trend towards lower values has been observed after storage for 12 h (0.79 [0.65–0.84], p = 0.112). Furthermore, cytochrome *c* control efficiency, which reflects the integrity of the mitochondrial outer membrane, revealed a slight increase at 8 h and 12 h HTS, indicating higher membrane permeability (p = 0.199 and p = 0.112) (**Figure 5**).

Real-Time Confocal Microscopy and Conventional Histology

In contrast to mitochondrial respiration, tissue integrity and cell viability could be preserved for up to 12 h of HTS, as demonstrated by RTCM analysis (**Figure 6**). The semiquantitative RTCM score revealed no changes during the storage period (data not shown). The microscopic evaluation of H&E sections also showed no differences regarding histopathology when comparing fresh biopsies with 12 h HTS biopsies (**Supplementary Figure S1**). This also underlines the sensitivity of HRR in comparison to other assessment methods during cold storage.

DISCUSSION

Herein we investigated the feasibility of prolonged hypothermic storage of tissue biopsies taken during NMP for later HRR

analysis. Testing of mitochondrial respiration in livers undergoing NMP can provide information on organ function and displays a reliable biomarker for decision-making [6, 17]. The HRR parameters cytochrome c control efficiency, *P-L* control efficiency, and LEAK respiration of the succinate-linked pathway correlated with the clinical outcome after transplantation in a previous clinical study. While HRR is a sensitive method, it requires technical proficiency of the operators and measurement immediately after sample collection. Enabling postponement of mitochondrial respiratory assessment would allow to analyze more biopsies and reduce the burden on personnel. This aspect is particularly relevant in scenarios where a higher number of biopsies need to be assessed for a comprehensive understanding of organ viability.

Our findings shed light on the feasibility of preservation of mitochondrial function in hypothermically stored liver biopsies. It indicates that prolonged periods of cold storage following liver NMP might be harmful, while shorter storage times of up 4 h may be acceptable. Our results demonstrated that mitochondrial respiration remained stable within this timeframe, and thus reliable results on mitochondrial function can be obtained after this storage period. This could guide decision-making processes pertaining to organ acceptance or utilization. Interestingly, beyond the 4-h mark mitochondrial respiration and efficiency tend to decrease, while the mitochondrial outer membrane damage increases. This contrasts earlier findings of our department [15], where the preservation of OXPHOS capacity and cytochrome c control efficiency during 48 h of hypothermic storage in porcine livers were analyzed. The first key difference between these experiments is that livers analyzed in the present study have variable and notably longer cold ischemic



FIGURE 5 | Mitochondrial respiratory states and coupling control efficiencies. LEAK respiration (A), OXPHOS capacity (B), OXPHOS(c) capacity (C), Electron transfer capacity (D), P-L control efficiency (E) to evaluate efficiency of ATP production, and cytochrome c control efficiency (F) to evaluate integrity of the outer mitochondrial membrane. Data are represented as median, interquartile range and min-max. Friedman's test with Dunn's multiple comparisons was applied to assess differences compared to 0 h samples. No significant differences were found ($\rho < 0.05$).

TABLE 4 Coupling control efficiencies of the succinate-linked pathway.							
Coupling control efficiencies	0 h	4 h	p-value	8 h	p-value	12 h	<i>p</i> -value
P-L control efficiency	0.85 [0.83–0.86]	0.81 [0.80–0.88]	>0.999	0.81 [0.78–0.83]	0.150	0.79 [0.65–0.84]	0.112
Cytochrome c control efficiency	0.07 [0.04–0.09]	0.03 [0.01–0.10]	>0.999	0.11 [0.08–0.18]	0.199	0.11 [0.08–0.26]	0.112

Data are represented as median, interquartile range and min-max. Friedman's test with Dunn's multiple comparisons was applied to assess differences compared to 0 h samples.

times and have undergone reperfusion in the form of NMP, while in the previous study no machine perfusion was applied. This finding is in line with a recent study of Horváth *et al.* [22]. They investigated the effect of prolonged CIT and subsequent normothermic isolated perfusion on mitochondrial function in DCD rat livers. Mitochondrial respiration was decreased after



after normothermic machine perfusion and after 12 h hypothermic storage. Images from real-time confocal microscopy portray cells stained with Hoechst 33,342 (in blue, detection of all nuclei), PI (in red, detection of dead cell's nuclei), WGA (in green, cell morphology), and the merged view. The scale bars denote 100 µm.

24 h of additional CIT compared to the control group, although in their setting no oxygen carrier was applied, which negatively affects mitochondrial metabolism under normothermic conditions [23]. Furthermore, it has already been shown in small animal models that NMP leads to altered metabolomic signatures in liver tissue which affect mitochondrial biochemical pathways, such as reduced succinate levels. This holds also potential for using different perfusion models and various optimizing supplements to positively influence liver metabolism during NMP [24, 25]. The second major difference is the quality of the organs. Healthy porcine livers without preexisting damage seem to be less sensitive to mitochondrial damage when compared to human livers coming from aged, ECD and/or DCD donors. While in our recent study no decrease in mitochondrial respiration during the first 6 h of NMP was found in the overall cohort, individual changes in some livers could be demonstrated.

While the decreased respiration in biopsies at both the 8 h and 12 h HTS preservation did not reach statistical significance, the trends are clear, and this nonetheless underscores a potential threat of prolonged HTS following NMP. The decrease in mitochondrial function after reaching 8 h of HTS preservation could be attributed to several underlying factors such as the accumulation of metabolic byproducts, disruption of cellular homeostasis and ion gradients, oxidative stress, and consequent alterations in mitochondrial membrane integrity [1, 11]. This level of degradation at 8 h may warrant further exploration to better understand the mechanisms involved and potentially improve preservation methods. It is noteworthy that both the RTCM and the histopathological examinations did not show any appreciable differences between samples at 0 and 12 h of HTS.

Our study acknowledges certain limitations that pertain to both HRR and the organ quality. The reproducibility of HRR measurements can be influenced by sample variability, experimental conditions, technical factors, and from the measurement process itself, often influenced by the individual conducting the measurements. This can impact accuracy of the results, highlighting the importance of standardized protocols, training, quality control, and minimizing inter-operator variability [26]. In addition, the relatively small sample number of N = 5 livers and singlecenter design, indicate a possible frailty, which needs to be proven by further centers. While the stability of mitochondrial respiration up to 4 h upon sample collection in an NMP setting could be demonstrated, the actual clinical outcome could not be predicted in the 4 h samples as 4 of 5 livers have been rejected during NMP, based on conventional organ quality criteria.

In summary, this study accentuates the value of HRR in assessing the bioenergetic efficacy of livers during cold storage and the dynamics of mitochondrial stability and function during static cold storage following NMP. HTS of tissue biopsies of livers undergoing NMP can be performed for no longer than 4 h. While further research is needed, our findings underscore the potential of this simple method as a tool for informed decision-making during evaluation for organ transplantation, with the added benefit of easing the logistic challenges of HRR and facilitating concurrent assessment of a larger number of biopsies.

DATA AVAILABILITY STATEMENT

The raw data supporting the conclusion of this article will be made available by the corresponding author upon request.

ETHICS STATEMENT

The studies involving humans were approved by the Ethics Commission of the Medical University of Innsbruck. The studies were conducted in accordance with the local legislation and institutional requirements. The participants provided their written informed consent to participate in this study.

AUTHOR CONTRIBUTIONS

JH: study design, data acquisition, data analysis, interpretation of the data, drafting; AK: study design, data acquisition, data analysis, interpretation of the data, drafting; MS: study design, data analysis, and revising the manuscript; MLB: data analysis, and revising the manuscript; MH: data analysis, interpretation of the data, and revising the manuscript; BZ: data analysis, interpretation of the data, and revising the manuscript; DÖ: funding acquisition, revising the manuscript; RO: revising the manuscript, analysis and interpretation of clinical data, TH: revising the manuscript; SS: funding acquisition, revising the manuscript; ATM: study design, data analysis, interpretation of the data, drafting and revising the manuscript. All authors contributed to the article and approved the submitted version.

REFERENCES

- Nasralla D, Coussios CC, Mergental H, Akhtar MZ, Butler AJ, Ceresa CDL, et al. A Randomized Trial of Normothermic Preservation in Liver Transplantation. *Nature* (2018) 557(7703):50–6. doi:10.1038/s41586-018-0047-9
- Verhoeven CJ, Farid WR, de Jonge J, Metselaar HJ, Kazemier G, van der Laan LJ. Biomarkers to Assess Graft Quality During Conventional and Machine Preservation in Liver Transplantation. J Hepatol (2014) 61(3):672–84. doi:10. 1016/j.jhep.2014.04.031
- van Leeuwen OB, de Vries Y, Fujiyoshi M, Nijsten MWN, Ubbink R, Pelgrim GJ, et al. Transplantation of High-Risk Donor Livers After *Ex Situ* Resuscitation and Assessment Using Combined Hypo⁻ and Normothermic Machine Perfusion: A Prospective Clinical Trial. *Ann Surg* (2019) 270(5): 906–14. doi:10.1097/SLA.000000000003540
- Weissenbacher A, Bogensperger C, Oberhuber R, Meszaros A, Gasteiger S, Ulmer H, et al. Perfusate Enzymes and Platelets Indicate Early Allograft Dysfunction After Transplantation of Normothermically Preserved Livers. *Transplantation* (2022) 106:792–805. doi:10.1097/TP.000000000003857
- Mergental H, Laing RW, Kirkham AJ, Perera M, Boteon YL, Attard J, et al. Transplantation of Discarded Livers Following Viability Testing With Normothermic Machine Perfusion. *Nat Commun* (2020) 11(1):2939. doi:10. 1038/s41467-020-16251-3
- Meszaros AT, Hofmann J, Buch ML, Cardini B, Dunzendorfer-Matt T, Nardin F, et al. Mitochondrial Respiration During Normothermic Liver Machine Perfusion Predicts Clinical Outcome. *eBioMedicine* (2022) 85:104311. doi:10. 1016/j.ebiom.2022.104311
- 7. Schlegel A, Muller X, Mueller M, Stepanova A, Kron P, de Rougemont O, et al. Hypothermic Oxygenated Perfusion Protects From Mitochondrial Injury

FUNDING

The author(s) declare that financial support was received for the research, authorship, and/or publication of this article. The work was funded by the "In Memoriam Dr. Gabriel Salzner Stiftung", awarded to SS and DÖ.

CONFLICT OF INTEREST

The authors declare that the research was conducted in the absence of any commercial or financial relationships that could be construed as a potential conflict of interest.

ACKNOWLEDGMENTS

The authors thank all members of the NMP team at the organLife Laboratory, Department of Visceral, Transplant and Thoracic Surgery. We thank Oroboros Instruments for providing technical support for high-resolution respirometry.

SUPPLEMENTARY MATERIAL

The Supplementary Material for this article can be found online at: https://www.frontierspartnerships.org/articles/10.3389/ti.2024. 12787/full#supplementary-material

Before Liver Transplantation. *EBioMedicine* (2020) 60:103014. doi:10.1016/ j.ebiom.2020.103014

- Darius T, Vergauwen M, Maistriaux L, Evrard R, Schlegel A, Mueller M, et al. Intermittent Surface Oxygenation Results in Similar Mitochondrial Protection and Maintenance of Aerobic Metabolism as Compared to Continuous Oxygenation During Hypothermic Machine Kidney Machine Perfusion. *J Clin Med* (2023) 12(11):3731. doi:10.3390/jcm12113731
- Panconesi R, Carvalho MF, Eden J, Fazi M, Ansari F, Mancina L, et al. Mitochondrial Injury During Normothermic Regional Perfusion (NRP) and Hypothermic Oxygenated Perfusion (HOPE) in a Rodent Model of DCD Liver Transplantation. *EBioMedicine* (2023) 98:104861. doi:10.1016/j.ebiom.2023. 104861
- Kesseli SJ, Gloria JN, Abraham N, Halpern SE, Cywinska GN, Zhang M, et al. Point-of-Care Assessment of DCD Livers During Normothermic Machine Perfusion in a Nonhuman Primate Model. *Hepatol Commun* (2021) 5(9): 1527–42. doi:10.1002/hep4.1734
- Horvath T, Jasz DK, Barath B, Poles MZ, Boros M, Hartmann P. Mitochondrial Consequences of Organ Preservation Techniques During Liver Transplantation. *Int J Mol Sci* (2021) 22(6):2816. doi:10.3390/ ijms22062816
- Garcia-Roche M, Casal A, Carriquiry M, Radi R, Quijano C, Cassina A. Respiratory Analysis of Coupled Mitochondria in Cryopreserved Liver Biopsies. *Redox Biol* (2018) 17:207–12. doi:10.1016/j.redox.2018.03.008
- Steinlechner-Maran R, Eberl T, Kunc M, Schröcksnadel H, Margreiter R, Gnaiger E. Respiratory Defect as an Early Event in Preservation-Reoxygenation Injury of Endothelial Cells. *Transplantation* (1997) 63(1): 136–42. doi:10.1097/00007890-199701150-00025
- Kuznetsov AV, Brandacher G, Steurer W, Margreiter R, Gnaiger E. Estimation of Mitochondrial Damage in Heart Preservation. *Transpl Proc* (1997) 31(1-2): 992. doi:10.1016/s0041-1345(98)01871-5

- Kuznetsov AV, Strobl D, Ruttmann E, Konigsrainer A, Margreiter R, Gnaiger E. Evaluation of Mitochondrial Respiratory Function in Small Biopsies of Liver. *Anal Biochem* (2002) 305(2):186–94. doi:10.1006/abio. 2002.5658
- Cardini B, Oberhuber R, Fodor M, Hautz T, Margreiter C, Resch T, et al. Clinical Implementation of Prolonged Liver Preservation and Monitoring Through Normothermic Machine Perfusion in Liver Transplantation. *Transplantation* (2020) 104(9):1917–28. doi:10.1097/ TP.000000000003296
- Hofmann J, Meszaros AT, Buch ML, Nardin F, Hackl V, Strolz CJ, et al. Bioenergetic and Cytokine Profiling May Help to Rescue More DCD Livers for Transplantation. *Int J Mol Sci* (2023) 24(11):9536. doi:10.3390/ijms24119536
- Doerrier C, Garcia-Souza LF, Krumschnabel G, Wohlfarter Y, Meszaros AT, Gnaiger E. High-Resolution FluoRespirometry and OXPHOS Protocols for Human Cells, Permeabilized Fibers From Small Biopsies of Muscle, and Isolated Mitochondria. *Methods Mol Biol* (2018) 1782:31–70. doi:10.1007/ 978-1-4939-7831-1_3
- Gnaiger E. Mitochondrial Pathways and Respiratory Control. An Introduction to OXPHOS Analysis. *Bioenergetic Commun* (2020) 2. doi:10.26124/bec:2020-0002
- Cardini B, Fodor M, Hermann M, Wieser V, Hautz T, Mellitzer V, et al. Live Confocal Imaging as a Novel Tool to Assess Liver Quality: Insights From a Murine Model. *Transplantation* (2020) 104(12):2528–37. doi:10.1097/TP. 000000000003405
- Weissenbacher A, Oberhuber R, Hermann M, Soleiman A, Ulmer H, Mayer G, et al. Live Confocal Tissue Assessment With SYTO16/PI and WGA Staining Visualizes Acute Organ Damage and Predicts Delayed Graft Function in Kidney Transplantation. *Ann Surg* (2019) 270(5):915–22. doi:10.1097/SLA. 000000000003511

- 22. Horvath T, Sandor L, Barath B, Donka T, Barath B, Mohacsi A, et al. Methane Admixture Protects Liver Mitochondria and Improves Graft Function After Static Cold Storage and Reperfusion. *Antioxidants (Basel)* (2023) 12(2):271. doi:10.3390/antiox12020271
- Dondossola D, Santini A, Lonati C, Zanella A, Merighi R, Vivona L, et al. Human Red Blood Cells as Oxygen Carriers to Improve *Ex-Situ* Liver Perfusion in a Rat Model. *J Clin Med* (2019) 8(11):1918. doi:10.3390/ jcm8111918
- 24. Lonati C, Dondossola D, Zizmare L, Battistin M, Wüst L, Vivona L, et al. Quantitative Metabolomics of Tissue, Perfusate, and Bile From Rat Livers Subjected to Normothermic Machine Perfusion. *Biomedicines* (2022) 10(3): 538. doi:10.3390/biomedicines10030538
- Guo Z, Zhan L, Gao N, Zhang Z, Huang S, Wang L, et al. Metabolomics Differences of the Donor Livers Between *in Situ* and *Ex Situ* Conditions During Ischemia-Free Liver Transplantation. *Transplantation* (2023) 107(5): e139–e151. doi:10.1097/TP.00000000004529
- Baglivo E, Cardoso LH, Cecatto C, Gnaiger E. Statistical Analysis of Instrumental Reproducibility as Internal Quality Control in High-Resolution Respirometry. *Bioenerg Commun* (2022) 8. doi:10.26124/bec: 2022-0008

Copyright © 2024 Hofmann, Kofler, Schartner, Buch, Hermann, Zelger, Öfner, Oberhuber, Hautz, Schneeberger and Meszaros. This is an open-access article distributed under the terms of the Creative Commons Attribution License (CC BY). The use, distribution or reproduction in other forums is permitted, provided the original author(s) and the copyright owner(s) are credited and that the original publication in this journal is cited, in accordance with accepted academic practice. No use, distribution or reproduction is permitted which does not comply with these terms.

GLOSSARY

ATP	adenosine triphosphate
СССР	Carbonyl cyanide p-trifluoromethoxyphenyl hydrazone
СІТ	cold ischemia time
CVA	cerebrovascular accident
DBD	donation after brain death
DCD	donation after circulatory death
DRI	donor risk index
ECD	extended criteria donor
ET	electron transfer
F	FAO-linked pathway
FAO	fatty acid oxidation
FCR	flux control ratio
FMN	flavin mononucleotide
HRR	high-resolution respirometry
HTS	hypothermic storage
IQR	interquartile range
L	LEAK respiration
LEAK	the state of non-phosphorylating resting respiration in the absence
	of ADP
LT	liver transplantation
N	NADH-linked pathway
NADH	nicotinamide adenine dinucleotide hydrogen
NMP	normothermic machine perfusion
OXPHOS	oxidative phosphorylation
Р	OXPHOS: Respiratory capacity of mitochondria in the ADP-activated
	state of oxidative phosphorylation, at saturating concentrations of ADP
	and inorganic phosphate, oxygen, and defined substrates
PI	propidium iodide
ROX	residual oxygen consumption
RTCM	real-time confocal microscopy
S	succinate-linked pathway
SUIT	Substrate-Uncoupler-Inhibitor Titration
WGA	wheat germ agglutinin